

Environmental conditions




Operating temperature	15-40°C (59-104°F)
Storage temperature	5-45°C (41-113°F)
Humidity	20-80%, non-condensing
Max pressure	7 bars (102 psi)

Fluidic characteristics

Tube port fittings	Standard ¼-28 UNF, flat-bottom
Wetted materials	PTFE, PCTFE and borosilicate glass
Channel diameter	0.5 mm (0.020 in) / 0.3-1.6 mm (0.012-0.063 in) available upon request
Axial port (stainless steel)	¼-28 UNF, flat-bottom (10 mm depth)

Mechanical characteristics

VALVE MODELS

Series	DISTRIBUTION	ON/OFF	SWITCH
Description	Choose 1 among N ports	Simple switch: flow or no flow	Channels are linked two by two
Liquid path			

For known aliases and catalogue models, please refer to <https://amf.ch/product/rvm-microfluidic-electric-rotary-valve/>

Glossary reminder

Internal volume:	Volume inside the system, from entrance to exit.
Dead volume:	Volume that is “stuck” in the system (dead end), which is not cleanly swept and relies on diffusion to clear out.
Carryover volume:	Volume of liquid that will be mixed with the next liquid. It is not stuck, but will be swept next time a liquid passes.

RVM Series - Industrial Microfluidic Rotary Valve

DATASHEET





MINI MODULE (Ref. P202-O)

Weight (without valve head)	61 g
Dimensions	Ø32 x 57 mm
External fixation system	3x M3 screws, position shown in orange
Rotation time for 180°	Down to 600 ms
Switching time Port to Port	Down to 180 ms



Fixations for MINI MODULE (options)



RVM Series - Industrial Microfluidic Rotary Valve


DATASHEET



Advanced Microfluidics

FAST MODULE (Ref. P201-O)

Weight (without valve head)	442 g
Dimensions	42.3 x 63 x 101 mm
External fixation system	6x M3 screws, position shown in orange
Rotation time for 180°	Down to 400 ms
Switching time Port to Port	Down to 136 ms



Advanced Microfluidics SA
 Chemin de la Dent d'Oche 1a
 CH-1024 Ecublens, Switzerland
 Phone: +41 (0)21 552 14 30
 Email: info@amf.ch
 Web: www.amf.ch


RVM Series - Industrial Microfluidic Rotary Valve

DATASHEET



LOW POWER MODULE (Ref. P200-O)

Weight (without valve head)	245 g
Dimensions	29 x 41 x 112 mm
External fixation system	6x M3 screws, position shown in orange
Rotation time for 180°	1.5 s



Electrical characteristics

MINI MODEL (Ref. P202-O)

Power	22-32 VDC, 0.5 A peak
Required cables for operation	USB C male to USB A or C male and custom cable with Pico-Lock™ 1mm pitch 4 pos. male

FAST MODEL (Ref. P201-O)

Power	18-24 VDC, 2 A peak
Required cables for operation	USB mini B male to USB A male or custom cable with PicoBlade™ 4 pos. male, and custom power cord

LOW POWER MODEL (Ref. P200-O)

Power	5-10 VDC, 0.5 A peak
Required cables for operation	USB mini B male to USB A male or custom cable with PicoBlade™ 6 pos. male

DO NOT PLUG THE RVM INTO A USB HUB. Please plug it into a main USB port directly for better performance.

Communication interface

MINI MODEL (Ref. P202-O)

Interface	USB C, 3x Pico-Lock™ 1mm pitch 4 pos.
Communication type	Serial communication on USB C I2C + Power on 1st Pico-Lock™ TTL (or other upon request) on 2nd Pico-Lock™ RS232 (RS485 or CAN upon request) + Power on 3rd Pico-Lock™

FAST MODEL (Ref. P201-O)

Interface	USB mini B, PicoBlade 4 pos., DB9 (other upon request)
Communication type	Serial communication on USB mini I2C on PicoBlade™ RS232 on DB9 (RS485 on DB9 upon request) (other upon request)

LOW POWER MODEL (Ref. P200-O)

Interface	USB mini B, PicoBlade 6 pos.
Communication type	Serial communication on USB mini I2C on PicoBlade™ (other upon request)

Through the USB, the valve controller is seen as a virtual serial port. There is an on-board USB to serial chip, so the drivers are automatically installed on recent Windows 7 and above.

SERIAL COMMUNICATION PARAMETERS

Baud rate	9600 (57600 for P202-O by default)
Data Bit	8
Parity	None
Stop Bit	1
Handshake	None
End line	<CR>

View the RVM operating manual for the extensive command set.

Revisions

Version	Author	Release Date	Modifications	Revised By	Approved By
05.03	NCR	22/05/2026	P202-O integration	CAL	ARD

The information and data contained in this documentation is subject to change without notice. The reproduction, distribution, and utilization of this document as well as the communication of its contents to others without explicit authorization is prohibited.

All rights reserved in the event of the grant of a patent, utility model or design.

The general terms and conditions of Advanced Microfluidics SA apply. Alternative agreements must be in written form.

Advanced Microfluidics SA shall not be liable for errors contained in this document or for incidental or consequential damages in connection with the furnishing, performance, or use of this material.