


# TECHNICAL NOTE

## TWO SAMPLES, TWO ANALYSERS

Advanced  
MicroFluidics  
Innovative  
automation solutions



### WHAT ARE WE DOING?

Using a 4-port switch valve (VS1-4), it is possible to send two different samples through 2 detectors (one at a time) with a single setup. Automating an analysis procedure then becomes an option!

### WHAT'S IN IT FOR YOU?

- Repeatability
- Automated system
- Flow control
- Minimal bench space

### CONTACT US FOR MORE INFORMATION

Advanced Microfluidics SA  
Ch. de la Dent d'Oche 1a  
CH-1024 Ecublens  
Switzerland

[www.amf.ch](http://www.amf.ch)  
[info@amf.ch](mailto:info@amf.ch)  
+41 21 552 14 30