




TECHNICAL NOTE

LIQUID RECIRCULATION - EXTERNAL PUMP



WHAT ARE WE DOING?

Using a 2-position 6-port switch valve (VS1-6), it is possible to create a system that recirculates the liquid between two vials, all the while maintaining the flow in the same direction in your flow cell.

WHAT'S IN IT FOR YOU?

- Liquid economy
- Shear stress in the flow cell always in the same direction
- Automated system
- Flow control

CONTACT US FOR MORE INFORMATION

Advanced Microfluidics SA
Ch. de la Dent d'Oche 1a
CH-1024 Ecublens
Switzerland

www.amf.ch
info@amf.ch
+41 21 552 14 30