


# NOW AVAILABLE

do more with a  
single machine



## 10 & 12 PORTS

distribution  
valves

THESE VALVES ARE COMPATIBLE WITH ALL OUR STANDARD  
MODULES:

RVM, SPM, LSPONE

### Valves Specifications

REF.	CONFIGURATION	WETTED MATERIALS	INTERNAL VOL.	CARRYOVER VOL.	CHANNEL $\phi$
VFD1-10-050	10 ports	PCTFE, PE*	4.5 $\mu$ L	2.8 $\mu$ L	0.5 mm
VFD1-10-100	10 ports	PCTFE, PE*	18.1 $\mu$ L	11 $\mu$ L	1 mm
VFD1-12-050	12 ports	PCTFE, PE*	4.5 $\mu$ L	2.8 $\mu$ L	0.5 mm

\*other materials available upon request

IF YOU WOULD LIKE MORE  
INFORMATION OR ANOTHER  
VARIATION, PLEASE GET IN  
TOUCH WITH US.

Advanced Microfluidics SA  
Ch. de la Dent d'Oche 1a  
CH-1024 Ecublens  
Switzerland

+41 21 552 14 30  
info@amf.ch  
www.amf.ch